

BPS 24V Power Converter for ResMed S9 Series

Battery Power Solutions' 24V Power Converter was specifically designed for use with ResMed S9 Series PAP devices and the C-100 Battery Pack. This single component connection provides a clean, efficient conversion from 24V to 12V direct current.

CONNECTION INSTRUCTIONS

- Connect the male end of the 24V Power Converter to the DC Out Jack of the Battery Pack.

NOTE: You must connect the 24V Power Converter to the Battery Pack before you connect it to the ResMed S9 Series device.

- If you are using the Battery Pack in the dual configuration, connect the DC Power Pigtail Cord to the DC Out Jack of both Battery Packs and connect the 24V Power Converter to the end of the DC Power Pigtail Cord.
- Connect the female end of the 24V Power Converter to the DC input port of the ResMed S9 Series device.
- Press the ON/OFF button on the PAP device and the device should start functioning.



RESMED DEVICE WARNING

If you are using the Battery Pack to power any device manufactured by ResMed **you must** use either a power converter or inverter to complete the connection between the device and the Battery Pack. A power converter or inverter must be used because ResMed devices operate from a 24V direct current and their plugs may run opposite polarity from other manufacturer devices. **Do not connect the Battery Pack directly to any ResMed device as doing so will damage the PAP device.**