

BPS 150W Sine Wave Power Inverter

The 150W modified sine wave power inverter can be used in conjunction with the Freedom™ or C-100 CPAP Battery to:

- Safely and efficiently charge the CPAP battery via the cigarette lighter/power socket of a vehicle.
- Power certain PAP device models including ResMed S8 series, Fisher & Paykel, Transcend, 3B Medical, and Apex Medical devices.
- Power other small electronic devices such as cell phones and tablets.

The 150W modified sine wave power inverter is not proprietary to Battery Power Solutions CPAP batteries and can be used to power devices in a vehicle or in conjunction with other types of external batteries.

SPECIFICATIONS

Input: 12V DC

Output: 110 - 120V AC

USB Output: 5V DC, 1 A/2.1 A

Rated Power: 150W (200W max)

Size: 3.5" x 2.25" x 1" (L x W x H) Cord Length: 31.5" Weight: 6.2 oz


CHARGING THE FREEDOM™ OR C-100 CPAP BATTERY

- Plug the AC power supply into the DC IN jack of the battery.

Dual Battery Connection: Plug the AC power supply into the DC power pigtail cord included in the CPAP battery kit and then plug the pigtail cord into the DC IN jacks of both batteries.

- Plug the AC power supply into the inverter.
- Plug the inverter into the cigarette lighter/power socket of the vehicle.
- When charging is complete disconnect the inverter from the vehicle and the battery.

SINGLE BATTERY CONNECTION EXAMPLE



POWERING OTHER DEVICES

- Plug the DC input cord included with the CPAP battery kit into the DC OUT jack of the battery.

Dual Battery Connection: Plug the DC input cord into the DC power pigtail cord included with the CPAP battery kit and then insert the pigtail cord into the DC OUT jacks of both batteries.

- Plug the inverter into the DC input cord.
- Plug the AC power supply or USB cord of your device into the inverter.

DUAL BATTERY CONNECTION EXAMPLE

